


ATS - Autotrans

*Closed Low-Vacuum-System
for postoperative autologous retransfusion*



120-000

ATS-Autotrans
Autotransfusionssystem

▶ **PATIENT SAFE**

▶ **EASY OF USE**

▶ **COST EFFECTIVE**

ATS - Autotrans



210-000 ATS-40 Transfusionset
with 40µ-Filter



120-500 Wastebag

- Low-Vacuum – reduced Hemolysis
- Easy of Use
- One closed system
- No device (pump) needed
- Autologous blood reduce usage of homologues red cells
- 1000 ml Retransfusion bag
- 2 One-way-valves
- 200 µm Filter integrated in the system
- Cost effective
- Latex free

CE 0482

List. Nr.	Description	units/carton
120-000	ATS-Autotrans Autotransfusion System with 3 Drainage connector and de-connectable Drainage-extension, 375ml bellow and 1000ml Retransfusion bag with One-way-valve and 200µ - Filter	10
210-000	ATS-40 Transfusion Set with 40 µm –Filter, 150 cm, Roller-clamps, Luer Lock	60
120-500	Drainage-Bag Wastebag 700ml	60

STRADIS

Stradis HandelsgmbH
Silberbachweg 12
A-2011 Sierndorf
Tel./Fax: +43 2267 43091
mail: office@stradis.at
web: www.stradis.at